

PACO43430

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Acts on tyrosine phosphorylated proteins, low-MW aryl phosphates and natural and synthetic acyl phosphates. Isoform 3 does not possess phosphatase activity.

Gene ID:

ACP1

Uniprot

P24666

Synonyms:

Low molecular weight phosphotyrosine protein phosphatase (LMW-PTP) (LMW-PTPase) (EC 3.1.3.48) (Adipocyte acid, phosphatase) (Low molecular weight cytosolic acid, phosphatase) (EC 3.1.3.2) (Red cell acid, phosphatase 1), ACP1

Immunogen:

Recombinant Human Low molecular weight phosphotyrosine protein phosphatase protein (1-158AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

N/A

N/A