

AEBP2 Antibody



PACO43436

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

DNA-binding transcriptional repressor. May interact with and stimulate the activity of the PRC2 complex, which methylates 'Lys-9' and 'Lys-27' residues of histone H3.

Gene ID:

AEBP2

Uniprot

Q6ZN18

Synonyms:

Zinc finger protein AEBP2 (Adipocyte enhancer-binding protein 2) (AE-binding protein 2), AEBP2


Immunogen:

Recombinant Human Zinc finger protein AEBP2 protein (248-517AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO43436 at dilution of 1:100.