

AIPL1 Antibody



PACO43439

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be important in protein trafficking and/or protein folding and stabilization.

Gene ID:

AIPL1

Uniprot

Q9NZN9

Synonyms:

Aryl-hydrocarbon-interacting protein-like 1, AIPL1, AIPL2


Immunogen:

Recombinant Human Aryl-hydrocarbon-interacting protein-like 1 protein (1-384AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO43439 at dilution of 1:100.