

# AIPL1 Antibody



PACO43440

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:2000,  
IHC:1:20-1:200

**Protein Background:**

May be important in protein trafficking and/or protein folding and stabilization.

**Gene ID:**

AIPL1

**Uniprot**

Q9NZN9

**Synonyms:**

Aryl-hydrocarbon-interacting protein-like 1, AIPL1, AIPL2

**Immunogen:**


Recombinant Human Aryl-hydrocarbon-interacting protein-like 1 protein (1-384AA).

**Storage:**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot. All lanes: Aryl-hydrocarbon-interacting protein-like 1 antibody at 5 $\mu$ g/ml + Rat brain tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 44, 37, 42, 41 kDa. Observed band size: 44 kDa.



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO43440 at dilution of 1:100.