

PACO43455

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

DIRAS3

Uniprot

O95661

Synonyms:

GTP-binding protein Di-Ras3 (Distinct subgroup of the Ras family member 3) (Rho-related GTP-binding protein Rho), DIRAS3, ARHI NOEY2 RHO1


Immunogen:

Recombinant Human GTP-binding protein Di-Ras3 protein (1-229AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO43455 at dilution of 1:100.