

PACO43470

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Enhances incorporation of BMP1 in the fibronectin matrix of connective tissues, and subsequent proteolytic activation of lysyl oxidase LOX. Induces cell attachment and spreading and plays a role in cell adhesion. May play a role in extracellular matrix mineralization.

**Gene ID:**

POSTN

**Uniprot**

Q15063

**Synonyms:**

Periostin (PN) (Osteoblast-specific factor 2) (OSF-2), POSTN, OSF2

**Immunogen:**


Recombinant Human Periostin protein (557-836AA).

**Storage:**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human breast cancer using PACO43470 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO43470 at dilution of 1:100.