

PPP1R1B Antibody



PACO43471

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Inhibitor of protein-phosphatase 1.

Gene ID:

PPP1R1B

Uniprot

Q9UD71

Synonyms:

Protein phosphatase 1 regulatory subunit 1B (DARPP-32) (Dopamine- and cAMP-regulated neuronal phosphoprotein), PPP1R1B, DARPP32


Immunogen:

Recombinant Human Protein phosphatase 1 regulatory subunit 1B protein (1-168AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO43471 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human epityphlon tissue using PACO43471 at dilution of 1:100.