

SGK3 Antibody



PACO43478

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

SGK3

Uniprot

Q96BR1

Synonyms:

Serine/threonine-protein kinase Sgk3 (EC 2.7.11.1) (Cytokine-independent survival kinase) (Serum/glucocorticoid-regulated kinase 3) (Serum/glucocorticoid-regulated kinase-like), SGK3, CISK SGKL


Immunogen:

Recombinant Human Serine/threonine-protein kinase Sgk3 protein (1-150AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human cervical cancer using PACO43478 at dilution of 1:100.