

SIRT6 Antibody



PACO43479

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

SIRT6

Uniprot

Q8N6T7

Synonyms:

NAD-dependent protein deacetylase sirtuin-6 (EC 3.5.1. -) (Regulatory protein SIR2 homolog 6) (SIR2-like protein 6), SIRT6, SIR2L6

Immunogen:

Recombinant Human NAD-dependent protein deacetylase sirtuin-6 protein (1-355AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

N/A

N/A