

FGFBP1 Antibody



PACO43506

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Acts as a carrier protein that release fibroblast-binding factors (FGFs) from the extracellular matrix (EM) storage and thus enhance the mitogenic activity of FGFs. Enhances FGF2 signaling during tissue repair, angiogenesis and in tumor growth.

Gene ID:

FGFBP1

Uniprot

Q14512

Synonyms:

Fibroblast growth factor-binding protein 1 (FGF-BP) (FGF-BP1) (FGF-binding protein 1) (FGFBP-1) (17 kDa heparin-binding growth factor-binding protein) (17 kDa HBGF-binding protein) (HBp17), FGFBP1, FGFBP HBP17


Immunogen:

Recombinant Human Fibroblast growth factor-binding protein 1 protein (24-234AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO43506 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO43506 at dilution of 1:100.