

MCPH1 Antibody



PACO43527

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Implicated in chromosome condensation and DNA damage induced cellular responses. May play a role in neurogenesis and regulation of the size of the cerebral cortex.

Gene ID:

MCPH1

Uniprot

Q8NEM0

Synonyms:

Microcephalin, MCPH1


Immunogen:

Recombinant Human Microcephalin protein (536-835AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO43527 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO43527 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).