

PACO43635

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

PA2 catalyzes the calcium-dependent hydrolysis of the 2-acyl groups in 3-sn-phosphoglycerides, this releases glycerophospholipids and arachidonic acid, that serve as the precursors of signal molecules.

Gene ID:

PLA2G1B

Uniprot

P04054

Synonyms:

Phospholipase A2 (EC 3.1.1.4) (Group IB phospholipase A2) (Phosphatidylcholine 2-acylhydrolase 1B), PLA2G1B, PLA2 PLA2A PPLA2

Immunogen:

Recombinant Human Phospholipase A2 protein (23-148AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

N/A

N/A