

THBS3 Antibody



PACO43645

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Adhesive glycoprotein that mediates cell-to-cell and cell-to-matrix interactions. Can bind to fibrinogen, fibronectin, laminin and type V collagen.

Gene ID:

THBS3

Uniprot

P49746

Synonyms:

Thrombospondin-3, THBS3, TSP3


Immunogen:

Recombinant Human Thrombospondin-3 protein (787-956AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO43645 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human epityphlon tissue using PACO43645 at dilution of 1:100.