

VENTX Antibody



PACO43646

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved in ventralization.

Gene ID:

VENTX

Uniprot

O95231

Synonyms:

Homeobox protein VENTX (VENT homeobox homolog) (VENT-like homeobox protein 2), VENTX, HPX42B VENTX2


Immunogen:

Recombinant Human Homeobox protein VENTX protein (1-258AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human colon tissue using PACO43646 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO43646 at dilution of 1:100.