

# CLN3 Antibody



PACO43658

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Involved in microtubule-dependent, anterograde transport of late endosomes and lysosomes.

**Gene ID:**

CLN3

**Uniprot**

Q13286

**Synonyms:**

Battenin (Batten disease protein) (Protein CLN3), CLN3, BTS

**Immunogen:**

Recombinant Human Battenin protein (1-280AA).

**Storage:**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---

N/A

N/A