

ALOXE3 Antibody



PACO43676

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

ALOXE3

Uniprot

Q9BYJ1

Synonyms:

Hydroperoxide isomerase ALOXE3 (EC 5.4.4.7) (Epidermis-type lipoxygenase 3) (Epidermal LOX-3) (e-LOX-3) (eLOX-3) (Hydroperoxy icosatetraenoate dehydratase) (EC 4.2.1.152), ALOXE3


Immunogen:

Recombinant Human Hydroperoxide isomerase ALOXE3 protein (1-280AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO43676 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human skin tissue using PACO43676 at dilution of 1:100.