

# RASGRP3 Antibody



PACO43826

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Guanine nucleotide exchange factor (GEF) for Ras and Rap1.

**Gene ID:**

RASGRP3

**Uniprot**

Q8IV61

**Synonyms:**

Ras guanyl-releasing protein 3 (Calcium and DAG-regulated guanine nucleotide exchange factor III) (Guanine nucleotide exchange factor for Rap1), RASGRP3, GRP3 KIAA0846

**Immunogen:**

Recombinant Human Ras guanyl-releasing protein 3 protein (541-690AA).

**Storage:**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---

N/A

N/A