

PACO43849

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Functions as a component of the PCAF complex. The PCAF complex is capable of efficiently acetylating histones in a nucleosomal context. The PCAF complex could be considered as the human version of the yeast SAGA complex.

Gene ID:

TAF5L

Uniprot

O75529

Synonyms:

TAF5-like RNA polymerase II p300/CBP-associated factor-associated factor 65 kDa subunit 5L (PCAF-associated factor 65 beta) (PAF65-beta), TAF5L, PAF65B


Immunogen:

Recombinant Human TAF5-like RNA polymerase II p300/CBP-associated factor-associated factor 65 kDa subunit 5L protein (1-325AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO43849 at dilution of 1:100.