

STMN4 Antibody



PACO43878

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Exhibits microtubule-destabilizing activity.

Gene ID:

STMN4

Uniprot

Q9H169

Synonyms:

Stathmin-4 (Stathmin-like protein B3) (RB3), STMN4


Immunogen:

Recombinant Human Stathmin-4 protein (1-216AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO43878 at dilution of 1:100.