

PACO43893

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Component of the asymmetric unit membrane (AUM); a highly specialized biomembrane elaborated by terminally differentiated urothelial cells. May play an important role in normal bladder epithelial physiology, possibly in regulating membrane permeability of superficial umbrella cells or in stabilizing the apical membrane through AUM/cytoskeletal interactions.

Gene ID:

UPK1B

Uniprot

O75841

Synonyms:

Uroplakin-1b (UP1b) (Tetraspanin-20) (Tspan-20) (Uroplakin Ib) (UPIb), UPK1B, TSPAN20

Immunogen:

Recombinant Human Uroplakin-1b protein (110-230AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

N/A

N/A