

TACSTD2 Antibody



PACO43895

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May function as a growth factor receptor.

Gene ID:

TACSTD2

Uniprot

P09758

Synonyms:

Tumor-associated calcium signal transducer 2 (Cell surface glycoprotein Trop-2) (Membrane component chromosome 1 surface marker 1) (Pancreatic carcinoma marker protein GA733-1), TACSTD2, GA733-1 M1S1 TROP2


Immunogen:

Recombinant Human Tumor-associated calcium signal transducer 2 protein (27-270AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO43895 at dilution of 1:100.