

PCDHA6 Antibody



PACO43899

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Potential calcium-dependent cell-adhesion protein. May be involved in the establishment and maintenance of specific neuronal connections in the brain.

Gene ID:

PCDHA6

Uniprot

Q9UN73

Synonyms:

Protocadherin alpha-6 (PCDH-alpha-6), PCDHA6, CNRS2


Immunogen:

Recombinant Human Protocadherin alpha-6 protein (1-230AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO43899 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO43899 at dilution of 1:100.