

# HYLS1 Antibody



PACO43962

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:2000,  
IHC:1:20-1:200

**Protein Background:**

HYLS1, also named as HLS, is widely expressed in fetal tissue. Mutation of HYLS1 will cause the hydrolethalus syndrome (HLS) which is an autosomal recessive lethal malformation syndrome characterized by multiple developmental defects of fetus.

**Gene ID:**

HYLS1

**Uniprot**

Q96M11

**Synonyms:**

Hydrolethalus syndrome protein 1, HYLS1, HLS

**Immunogen:**


Recombinant Human Hydrolethalus syndrome protein 1 protein (70-299AA).

**Storage:**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot. All lanes: HYL51 antibody at 8 $\mu$ g/ml + HeLa whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 34 kDa. Observed band size: 34 kDa.



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO43962 at dilution of 1:100.