

GLYATL1 Antibody



PACO43964

Product Information

Size:

50ul

Reactivity:

Human, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:200-1:1000,
IHC:1:20-1:200

Protein Background:

Acyltransferase which transfers an acyl group to the N-terminus of glutamine. Can use phenylacetyl-CoA as an acyl donor.

Gene ID:

GLYATL1

Uniprot

Q969I3

Synonyms:

Glycine N-acyltransferase-like protein 1 (EC 2.3.1.68) (Acyl-CoA: glycine N-acyltransferase-like protein 1) (Glutamine N-acyltransferase), GLYATL1, GNAT


Immunogen:

Recombinant Human Glycine N-acyltransferase-like protein 1 protein (1-302AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: GLYATL1 antibody at 8 μ g/ml + Rat liver tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 36, 39 kDa. Observed band size: 36 kDa.



Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO43964 at dilution of 1:100.