

APOBEC3C Antibody



PACO44071

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

APOBEC3C

Uniprot

Q9NRW3

Synonyms:

DNA dC->dU-editing enzyme APOBEC-3C (A3C) (EC 3.5.4. -) (APOBEC1-like) (Phorbolins I), APOBEC3C, APOBEC1L PBI


Immunogen:

Recombinant Human DNA dC->dU-editing enzyme APOBEC-3C protein (1-190AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44071 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human liver tissue using PACO44071 at dilution of 1:100.