

DAPK1 Antibody



PACO44116

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

DAPK1

Uniprot

P53355

Synonyms:

Death-associated protein kinase 1 (DAP kinase 1) (EC 2.7.11.1), DAPK1, DAPK

Immunogen:

Recombinant Human Death-associated protein kinase 1 protein (1141-1430AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

N/A

N/A