

HOXA1 Antibody



PACO44134

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. Acts on the anterior body structures. Seems to act in the maintenance and/or generation of hindbrain segments.

Gene ID:

HOXA1

Uniprot

P49639

Synonyms:

Homeobox protein Hox-A1 (Homeobox protein Hox-1F), HOXA1, HOX1F


Immunogen:

Recombinant Human Homeobox protein Hox-A1 protein (75-205AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO44134 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon tissue using PACO44134 at dilution of 1:100.