

HYAL1 Antibody



PACO44135

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:

May have a role in promoting tumor progression. May block the TGFB1-enhanced cell growth.

Gene ID:

HYAL1

Uniprot

Q12794

Synonyms:

Hyaluronidase-1 (Hyal-1) (EC 3.2.1.35) (Hyaluronoglucosaminidase-1) (Lung carcinoma protein 1) (LuCa-1), HYAL1, LUCA1

Immunogen:

Recombinant Human Hyaluronidase-1 protein (136-435AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: HYAL1 antibody at 2 μ g/ml + HepG2 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 49, 46, 28, 24, 20, 11, 38 kDa. Observed band size: 49 kDa.



Immunohistochemistry of paraffin-embedded human liver tissue using PACO44135 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO44135 at dilution of 1:100.