

SIVA1 Antibody



PACO44245

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Induces CD27-mediated apoptosis. Inhibits BCL2L1 isoform Bcl-x(L) anti-apoptotic activity. Inhibits activation of NF-kappa-B and promotes T-cell receptor-mediated apoptosis.

Gene ID:

SIVA1

Uniprot

O15304

Synonyms:

Apoptosis regulatory protein Siva (CD27-binding protein) (CD27BP), SIVA1, SIVA


Immunogen:

Recombinant Human Apoptosis regulatory protein Siva protein (1-110AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO44245 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44245 at dilution of 1:100.