

PACO44307

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Muscle assembly regulating factor. Mediates the antiparallel assembly of titin (TTN) molecules at the sarcomeric Z-disk.

**Gene ID:**

TCAP

**Uniprot**

O15273

**Synonyms:**

Telethonin (Titin cap protein), TCAP

**Immunogen:**


Recombinant Human Telethonin protein (1-167AA).

**Storage:**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO44307 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human heart tissue using PACO44307 at dilution of 1:100.