

PACO44375

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Involved in the suppression of bile acid, biosynthesis through down-regulation of CYP7A1 expression, following positive regulation of the JNK and ERK1/2 cascades. Stimulates glucose uptake in adipocytes. Activity requires the presence of KLB and FGFR4.

Gene ID:

FGF19

Uniprot

O95750

Synonyms:

Fibroblast growth factor 19 (FGF-19), FGF19


Immunogen:

Recombinant Human Fibroblast growth factor 19 protein (25-216AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO44375 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO44375 at dilution of 1:100.