

# PDE6B Antibody



PACO44405

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

This protein participates in processes of transmission and amplification of the visual signal. Necessary for the formation of a functional phosphodiesterase holoenzyme.

**Gene ID:**

PDE6B

**Uniprot**

P35913

**Synonyms:**

Rod cGMP-specific 3',5'-cyclic phosphodiesterase subunit beta (GMP-PDE beta) (EC 3.1.4.35), PDE6B, PDEB

**Immunogen:**

Recombinant Human Rod cGMP-specific 3',5'-cyclic phosphodiesterase subunit &beta; protein (1-190AA).

**Storage:**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---

N/A

N/A