

CCNE2 Antibody



PACO44437

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Essential for the control of the cell cycle at the late G1 and early S phase.

Gene ID:

CCNE2

Uniprot

O96020

Synonyms:

G1/S-specific cyclin-E2, CCNE2


Immunogen:

Recombinant Human G1/S-specific cyclin-E2 protein (255-404AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44437 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human thyroid tissue using PACO44437 at dilution of 1:100.