

PACO44451

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Cell adhesion molecule contributing to the interactive process required for myelination in the central nervous system.

Gene ID:

OMG

Uniprot

P23515

Synonyms:

Oligodendrocyte-myelin glycoprotein, OMG, OMGP


Immunogen:

Recombinant Human Oligodendrocyte-myelin glycoprotein protein (190-420AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO44451 at dilution of 1:100.