

PROK2 Antibody



PACO44452

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May function as an output molecule from the suprachiasmatic nucleus (SCN) that transmits behavioral circadian rhythm. May also function locally within the SCN to synchronize output. Potently contracts gastrointestinal (GI) smooth muscle.

Gene ID:

PROK2

Uniprot

Q9HC23

Synonyms:

Prokineticin-2 (PK2) (Protein Bv8 homolog), PROK2, BV8


Immunogen:

Recombinant Human Prokineticin-2 protein (20-129AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44452 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO44452 at dilution of 1:100.