

PACO44475

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Supplies the nucleotide substrate for thymidylate synthetase.

Gene ID:

DCTD

Uniprot

P32321

Synonyms:

Deoxycytidylate deaminase (EC 3.5.4.12) (dCMP deaminase), DCTD


Immunogen:

Recombinant Human Deoxycytidylate deaminase protein (1-178AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO44475 at dilution of 1:100.