

# ABCA3 Antibody



PACO44483

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Plays an important role in the formation of pulmonary surfactant, probably by transporting lipids such as cholesterol.

**Gene ID:**

ABCA3

**Uniprot**

Q99758

**Synonyms:**

ATP-binding cassette sub-family A member 3 (ABC-C transporter) (ATP-binding cassette transporter 3) (ATP-binding cassette 3), ABCA3, ABC3

**Immunogen:**


Recombinant Human ATP-binding cassette sub-family A member 3 protein (43-203AA).

**Storage:**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human lung tissue using PACO44483 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44483 at dilution of 1:100.