

TOB2 Antibody



PACO44554

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Anti-proliferative protein inhibits cell cycle progression from the G0/G1 to S phases.

Gene ID:

TOB2

Uniprot

Q14106

Synonyms:

Protein Tob2 (Protein Tob4) (Transducer of erbB-2 2), TOB2, KIAA1663 TOB4 TROB2


Immunogen:

Recombinant Human Protein Tob2 protein (1-70AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO44554 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO44554 at dilution of 1:100.