

TRIM69 Antibody



PACO44559

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May have E3 ubiquitin-protein ligase activity. May play a role in apoptosis.

Gene ID:

TRIM69

Uniprot

Q86WT6

Synonyms:

E3 ubiquitin-protein ligase TRIM69 (EC 2.3.2.27) (RFP-like domain-containing protein trimless) (RING finger protein 36) (RING-type E3 ubiquitin transferase TRIM69) (Tripartite motif-containing protein 69), TRIM69, RNF36


Immunogen:

Recombinant Human E3 ubiquitin-protein ligase TRIM69 protein (1-341AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO44559 at dilution of 1:100.