

ELF5 Antibody



PACO44573

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Transcriptionally activator that may play a role in regulating the later stages of keratinocytes terminal differentiation.

Gene ID:

ELF5

Uniprot

Q9UKW6

Synonyms:

ETS-related transcription factor Elf-5 (E74-like factor 5) (Epithelium-restricted ESE-1-related Ets factor) (Epithelium-specific Ets transcription factor 2) (ESE-2), ELF5, ESE2


Immunogen:

Recombinant Human ETS-related transcription factor Elf-5 protein (1-170AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO44573 at dilution of 1:100.