

MYO5C Antibody



PACO44661

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved in transferrin trafficking. Likely to power actin-based membrane trafficking in many physiologically crucial tissues.

Gene ID:

MYO5C

Uniprot

Q9NQX4

Synonyms:

Unconventional myosin-Vc, MYO5C


Immunogen:

Recombinant Human Unconventional myosin-Vc protein (270-420AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO44661 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO44661 at dilution of 1:100.