

NBR1 Antibody



PACO44663

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Acts probably as a receptor for selective autophagosomal degradation of ubiquitinated targets.

Gene ID:

NBR1

Uniprot

Q14596

Synonyms:

Next to BRCA1 gene 1 protein (Cell migration-inducing gene 19 protein) (Membrane component chromosome 17 surface marker 2) (Neighbor of BRCA1 gene 1 protein) (Protein 1A1-3B), NBR1, 1A13B KIAA0049 M17S2


Immunogen:

Recombinant Human Next to BRCA1 gene 1 protein (1-280AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO44663 at dilution of 1:100.