

PACO44691

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

ATP binding, magnesium ion binding, protein serine/threonine kinase activity, protein phosphorylation.

Gene ID:

SIK3

Uniprot

Q9Y2K2

Synonyms:

Serine/threonine-protein kinase SIK3 (EC 2.7.11.1) (Salt-inducible kinase 3) (SIK-3) (Serine/threonine-protein kinase QSK), SIK3, KIAA0999 QSK


Immunogen:

Recombinant Human Serine/threonine-protein kinase SIK3 protein (60-310AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human lung cancer using PACO44691 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human thymus tissue using PACO44691 at dilution of 1:100.