

SNX7 Antibody



PACO44693

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved in several stages of intracellular trafficking.

Gene ID:

SNX7

Uniprot

Q9UNH6

Synonyms:

Sorting nexin-7, SNX7


Immunogen:

Recombinant Human Sorting nexin-7 protein (85-330AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO44693 at dilution of 1:100.