

TAF11 Antibody



PACO44697

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Core TAFII present in both of the previously described TFIID species which either lack or contain TAFII30 (TFIID alpha and TFIID beta respectively).

Gene ID:

TAF11

Uniprot

Q15544

Synonyms:

Transcription initiation factor TFIID subunit 11 (TFIID subunit p30-beta) (Transcription initiation factor TFIID 28 kDa subunit) (TAF(II)28) (TAFII-28) (TAFII28), TAF11, TAF2I


Immunogen:

Recombinant Human Transcription initiation factor TFIID subunit 11 protein (1-211AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44697 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO44697 at dilution of 1:100.