

PACO44699

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Part of a specialized transcription system that mediates the transcription of most ribosomal proteins through the 5'-TCT-3' motif which is a core promoter element at these genes. Seems to also mediate the transcription of NF1. Does not bind the TATA box.

Gene ID:

TBPL1

Uniprot

P62380

Synonyms:

TATA box-binding protein-like protein 1 (TBP-like protein 1) (21 kDa TBP-like protein) (Second TBP of unique DNA protein) (STUD) (TATA box-binding protein-related factor 2) (TBP-related factor 2) (TBP-like factor) (TBP-related protein), TBPL1, TLF TLP TLP21 TRF2 TRP


Immunogen:

Recombinant Human TATA box-binding protein-like protein 1 protein (1-186AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44699 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human heart tissue using PACO44699 at dilution of 1:100.