

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

microtubule cytoskeleton organization, regulation of JNK cascade, regulation of MAPK cascade, regulation of neuron migration, regulation of neuron projection development, regulation of p38MAPK cascade, regulation of protein kinase C signaling.

Gene ID:

ULK4

Uniprot

Q96C45

Synonyms:

Serine/threonine-protein kinase ULK4 (EC 2.7.11.1) (Unc-51-like kinase 4), ULK4


Immunogen:

Recombinant Human Serine/threonine-protein kinase ULK4 protein (1-290AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44715 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human skin tissue using PACO44715 at dilution of 1:100.