

PACO44717

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

WDR47

Uniprot

O94967

Synonyms:

WD repeat-containing protein 47 (Neuronal enriched MAP-interacting protein) (Nemitin), WDR47, KIAA0893


Immunogen:

Recombinant Human WD repeat-containing protein 47 protein (350-600AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44717 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human thyroid tissue using PACO44717 at dilution of 1:100.