

ZBTB48 Antibody



PACO44726

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Binds to and regulates the J and/or S elements in MHC II promoter.

Gene ID:

ZBTB48

Uniprot

P10074

Synonyms:

Telomere zinc finger-associated protein (TZAP) (Krueppel-related zinc finger protein 3) (hKR3) (Zinc finger and BTB domain-containing protein 48) (Zinc finger protein 855), ZBTB48, HKR3 TZAP ZNF855


Immunogen:

Recombinant Human Telomere zinc finger-associated protein (1-280AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human spleen tissue using PACO44726 at dilution of 1:100.