

ZFAND3 Antibody




PACO44728


Product Information

Size:	Protein Background:
50ul	
Reactivity:	Gene ID:
Human	ZFAND3
Source:	Uniprot
Rabbit	Q9H8U3
Isotype:	Synonyms:
IgG	AN1-type zinc finger protein 3 (Testis-expressed protein 27), ZFAND3, TEX27
Applications:	Immunogen:
ELISA, IHC	Recombinant Human AN1-type zinc finger protein 3 protein (1-227AA).
Recommended dilutions:	Storage:
ELISA:1:2000-1:10000, IHC:1:20-1:200	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO44728 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO44728 at dilution of 1:100.